

INFORMACIÓN GENERAL

Nombre del producto : Embudo, acero inoxidable**Descripción :** Fabricado en acero inoxidable 316

DATOS TÉCNICOS

referencia	Øsuperior(mm)	Øvástago (mm)	h (mm)	unidades por ref.
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810



GENERAL INFORMATION

Product name : Funnel, stainless steel**Description :** Made of 316 stainless steel

TECHNICAL DATA

reference	Øtop (mm)	Østem (mm)	h (mm)	pcs/pack
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810

PRODUCT PHOTO



INFORMATIONS GÉNÉRALES

Nom produit : Entonnoir en acier inoxydable**Description :**

Fabriqué en acier inoxydable 316

DONNÉES TECHNIQUES

référence	Øsupérieur (mm)	Øtige (mm)	h (mm)	unités par ref.
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810



INFORMAZIONE GENERALE

Nome del prodotto : Imbuto in acciaio inossidabile**Descrizione :**

Realizzato in acciaio inossidabile 316

DATI TECNICI

referenza	Øsuperiore (mm)	Østantuffo (mm)	h (mm)	unità per ref.
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810



ALGEMENE INFORMATIE

Produktnaam : Roestvrijstalen trechter**Beschrijving :**

Gemaakt van 316 roestvrij staal

TECHNISCHE GEGEVENS

referentie	bovenØ(mm)	Østeel (mm)	h (mm)	stuks per ref.
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

VERPAKKING EN LOGISTIEKE GEGEVENS

Referentie	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810



ALLGEMEINE INFORMATIONEN

Produktname : Trichter aus Edelstahl
Beschreibung : Aus Edelstahl 316.

TECHNISCHE DATEN

Artikelnummer	Øoben(mm)	ØStiel(mm)	H (mm)	Stückzahlpro Artikel
FUNA-050-001	50	12	65	1
FUNA-060-001	60	15	65	1
FUNA-080-001	80	12	105	1
FUNA-100-001	100	15	120	1
FUNA-120-001	120	12	140	1
FUNA-150-001	150	20	170	1
FUNA-200-001	200	20	220	1
FUNA-250-001	250	19	230	1

VERPACKUNG UND LOGISTIKDATEN

Referenz	vol (l)	kg	TARIC	GTIN
FUNA-050-001	0,5	0,02	73269092	08434868109780
FUNA-060-001	0,252	0,02	73269092	08434868097650
FUNA-080-001	1	0,04	73269092	08434868109797
FUNA-100-001	1,2	0,06	73269092	08434868097667
FUNA-120-001	2	0,12	73269092	08434868109803
FUNA-150-001	3,825	0,18	73269092	08434868097674
FUNA-200-001	8,8	0,86	73269092	08434868097681
FUNA-250-001	10	0,94	73269092	08434868109810

