


0.5 ml

REFERENCE: PCRP-005-1K0



NOTES

1. All dimensions are in millimetres (mm).
2. Tolerances are $\pm$ as shown.
3. For general laboratory use.

SPECIFICATIONS


Reference:	PCRP-005-1K0
Volume:	0.5 ml
Material:	PP
Colour:	Transparent
Autoclavable:	Yes
Units per ref.:	1,000
Height:	$H = 31.2 \pm 0.5 \text{ mm}$
OD of mouth:	$7.87 \pm 0.3 \text{ mm}$

FEATURES

- Precise seal
- Accurate graduations
- Transparent material
- Special safe-lock lid

1.5 ml

REFERENCE: **PCRP-015-500**



NOTES

1. All dimensions are in millimetres (mm).
2. Tolerances are $\pm$ as shown.
3. For general laboratory use.


SPECIFICATIONS

Reference:	PCRP-015-500
Volume:	1.5 ml
Material:	PP
Colour:	Transparent
Autoclavable:	Yes
Units per ref.:	500
Height:	$H = 38.7 \pm 0.5$ mm
OD of mouth:	10.6 ± 0.3 mm

FEATURES

- Precise seal
- Accurate graduations
- Transparent material
- Special safe-lock lid


REFERENCES: PCRP-Y15-500 / PCRP-P15-500 / PCRP-B15-500 / PCRP-G15-500



NOTES

1. All dimensions are in millimetres (mm).
2. Tolerances are $\pm$ as shown.
3. Same dimensions apply to all colour variants.

COLOUR VARIANTS




SPECIFICATIONS

Volume:	1.5 ml
Material:	PP
Colour variants:	Yellow / Pink / Blue / Green
Autoclavable:	Yes
Units per ref.:	500
Height:	$H = 38.7 \pm 0.5 \text{ mm}$
OD of mouth:	$10.6 \pm 0.3 \text{ mm}$

FEATURES

- Precise seal
- Accurate graduations
- Material: PP
- Available in yellow, pink, blue and green
- Special safe-lock lid

REFERENCE: **PCRP-C20-500**



NOTES

1. All dimensions are in millimetres (mm).
2. Tolerances are $\pm$ as shown.
3. For general laboratory use.


SPECIFICATIONS

Reference:	PCRP-C20-500
Volume:	2.0 ml
Material:	PP
Colour:	Transparent
Autoclavable:	Yes
Units per ref.:	500
Height:	H = $42.15 \pm 0.5 \text{ mm}$
OD of mouth:	$12.8 \pm 0.2 \text{ mm}$
OD of tube:	$10.8 \pm 0.2 \text{ mm}$

FEATURES

- Precise seal
- Accurate graduations
- Transparent material
- Special safe-lock lid

REFERENCE: PCRP-050-200



NOTES

1. All dimensions are in millimetres (mm).
2. Tolerances are $\pm$ as shown.
3. For general laboratory use.

SPECIFICATIONS

Reference:	PCRP-050-200
Volume:	5.0 ml
Material:	PP
Colour:	Transparent
Autoclavable:	Yes
Units per ref.:	200
Height:	$H = 59.4 \pm 0.5$ mm
OD of mouth:	19.2 ± 0.3 mm
OD of tube:	16.6 ± 0.3 mm

FEATURES

- Precise seal
- Accurate graduations
- Transparent material
- Special safe-lock lid