

Product Engg. Drawing

DESSICATOR PLATE

DESSICATOR PLATE SPECIFICATIONS

Labbox Code	Outer Dia (D1) (mm)	Thickness (T1) (mm)	Center Hole Dia (D2) (mm)	Surface Holes Dia (mm)	No. Of Holes	Perforated Dia (mm) (D4)
DEPK-140-001	140	8	23	5	54	118
DEPK-190-001	190	8	23	5	120	167
DEPK-210-001	210	8	23	5	132	190
DEPK-240-001	240	8	26	5	152	220
DEPK-275-001	275	8	26	5	162	250
	± 2.5%	± 1mm	± 2 mm	± 0.25 mm		± 2.5 %

